



Performance Qualification

• *Injector and Flow Precision*

• *Instruments and Fluidics*

<i>Instrument Name</i>	<i>Model</i>	<i>Supplier's Name</i>	<i>Serial Number</i>
Pump	LPG-3400RS	Thermo Scientific	8043253
Autosampler	WPS-3000TSL	Thermo Scientific	8013202
Column Oven	TCC-3000SD	Thermo Scientific	Demo
UV Detector	DAD-3000RS	Thermo Scientific	8043573
Chromeleon Datasystem	V. 6.80 SR16 Build 5387	Thermo Scientific	40640

<i>Accessories</i>	<i>Description</i>
Back Pressure Device	Capillary (L:15 m; ID:0,18 mm)
Sample 4	Caffeine in Water 140 µg/ml
Solvent A	Water (HPLC-Grade)

• *Additional Information*

Customer: Customer's Name
 Operator: Operator's Name
 Operator's Jobtitle
 Execution Date: 2020-avr-09
 Next Qualification: 2020-oct

• *Limits, Values and Test Results*

	<i>Limit</i>	<i>Observed Value</i>	<i>Result</i>
Injector Precision - Area	<= 0,500 % RSD	0,419 % RSD	Test passed
Flow Precision - Ret. Time	<= 0,100 % RSD	0,000 % RSD	Test passed
	OR <= 0,0200 min SD	0,0000 min SD	

• *Signatures*

	<i>Signature</i>	<i>Date</i>
Submitter / Operator:	{seq.submitOperator.userName}	{seq.submitTime}
Reviewer:	{seq.reviewOperator.userName}	{seq.reviewTime}
Approver (e-sig. only):	{seq.approveOperator.userName}	{seq.approveTime}

• **Data for Injector and Flow Precision Test: Volume** **5,0 µl**

Sample Name	Ret.Time min Caffeine UV_VIS_1	Area mAU*min Caffeine UV_VIS_1
Injector and flow reproducibility_1	1,9950	115,036
Injector and flow reproducibility_2	1,9950	114,837
Injector and flow reproducibility_3	1,9950	115,070
Injector and flow reproducibility_4	1,9950	114,075
Injector and flow reproducibility_5	1,9950	114,034
Injector and flow reproducibility_6	1,9950	113,904
Injector and flow reproducibility_7	1,9950	113,923
Injector and flow reproducibility_8	1,9950	113,995
Injector and flow reproducibility_9	1,9950	114,459
Injector and flow reproducibility_10	1,9950	113,971
Average:	1,9950	114,330
RSD:	0,000 %	0,419 %
RSD Limit:	<= 0,100 %	<= 0,500 %
SD:	0,0000	
SD Limit:	<= 0,0200	
Result:	ok	ok

• ***Charts for Injector and Flow Precision Test***

• Chromatogram Overlay



Injection and Retention Time Precision

Name	Area	Ret.Time	Sample No.	Area
------	------	----------	------------	------

	mAU*min Caffeine UV_VIS_1	min Caffeine UV_VIS_1		Dev. from mean [%]
Injector and flow reproducibility_1	115,036	1,995	1	0,617
Injector and flow reproducibility_2	114,837	1,995	2	0,443
Injector and flow reproducibility_3	115,070	1,995	3	0,647
Injector and flow reproducibility_4	114,075	1,995	4	-0,224
Injector and flow reproducibility_5	114,034	1,995	5	-0,259
Injector and flow reproducibility_6	113,904	1,995	6	-0,373
Injector and flow reproducibility_7	113,923	1,995	7	-0,357
Injector and flow reproducibility_8	113,995	1,995	8	-0,294
Injector and flow reproducibility_9	114,459	1,995	9	0,113
Injector and flow reproducibility_10	113,971	1,995	10	-0,314
Average:	114,330	1,995		
RSD:	0,419 %	0,000 %		
SD:		0,000 min		

Calculation of RSD and Average

	Area:	Retention Time:
Sum of x:	1143,302879	19,95
sum of sq(x):	130716,2089	39,80025
Sq of sum (x):	1307141,474	398,0025
STD:	0,478604359	0
Average:	114,3302879	1,995
RSD:	0,004186155	0
Number of samples:	10	10

Determination of Pump Unit for Dionex DGPs


Sequence name: PQ_INJECTOR_FLOW_REPRO
Right end of the sequence name : EPRO
Pump's model number: LPG-3400RS
Pump's model variant: LPG
Pump unit:

Evaluation of Accessories (Following autosampler uses standard 3 or 5 instead of st

ModelNo	Loop Volume	Evaluation*)	Reference value for report
WPS-3000	40	0	S4
ACC-3000	200	0	
ACC-3000T	200	0	
WPS-3000	344	0	
WPS-3000	250	0	
Gina 50		0	
Gina 160		0	
Accela Sampler		0	
Accela Open Sampler		0	
OAS-3000		0	
VH-A10-A		0	
VR-A12-A		0	
VR-A13-A		0	
VR-A14-A		0	

VR-A15-A	0
VF-A10-A	0

*) ModelNos without loop volume setting uses different evaluation formulas compared to ModelNos with loop volume setting.







--

No.	Ret. Time
-----	-----------

	<i>Deviation from the mean value [%]</i>
1	0,000
2	0,000
3	0,000
4	0,000
5	0,000
6	0,000
7	0,000
8	0,000
9	0,000
10	0,000

Standard 4):