



Performance Qualification

• Step Accuracy: STD Gradient_3

• Instruments and Fluidics

Instrument Name	Model	Supplier's Name	Serial Number
Pump	LPG-3400RS	Thermo Scientific	8043253
UV Detector	DAD-3000RS	Thermo Scientific	8043573
Chromeleon Datasystem	V. 6.80 SR16 Build 5387	Thermo Scientific	40640

Accessories	Description
Back Pressure Device	Capillary (L:15 m; ID:0,18 mm)
Solvent A	Water (HPLC-Grade)
Solvent B for Gradient	Water + 0.1% Acetone

• Additional Information

Customer: Customer's Name
 Operator: Operator's Name
 Operator's Jobtitle
 Execution Date: 2020-avr-15
 Next Qualification: 2020-oct


• Limits and Test Results

	Limit	Observed max. Deviation	Result of all Steps
Step Accuracy	<= 2,000 %	0,771 %	Test passed
Step Ripple	<= 0,500 %	0,043 %	Test passed

• Signatures

	Signature	Date
Submitter / Operator:	{seq.submitOperator.userName}	{seq.submitTime}
Reviewer:	{seq.reviewOperator.userName}	{seq.reviewTime}
Approver (e-sig. only):	{seq.approveOperator.userName}	{seq.approveTime}

• **Chromatogram of STD Gradient_3**



Flow [ml/min]: 2,000

• **Data of STD Gradient_3**

<i>Observed Value [mAU]</i>	<i>Expected Value [%]</i>	<i>Calculated Value [%]</i>	<i>Abs. Critical Deviation [%]</i>	<i>Calculated Deviation [%]</i>	<i>Result</i>
0,00	0,00	0,000	2,000	0,000	ok
3,19	1,00	1,476	2,000	0,476	ok
109,65	50,00	50,771	2,000	0,771	ok
214,66	99,00	99,395	2,000	0,395	ok
215,97	100,00	100,000	2,000	0,000	ok

• **Ripple of STD Gradient_3**

<i>Step [%]</i>	<i>Ripple [mAU]</i>	<i>Calculated Ripple [%]</i>	<i>Critical Ripple [%]</i>	<i>Result</i>
1,00	0,094	0,043	0,500	ok
50,00	0,075	0,035	0,500	ok
99,00	0,067	0,031	0,500	ok

Definition:

Sample Name: STD Gradient_3

Gradienttype: STD

Chromeleon (c) 2017 Thermo Fisher Scientific Inc.
Version 6.80 SR16 Build 5387 (291597)

1

0

0

1

Gradient_accuracy_8043253 / PUMP_GRADIENT

Printed: 01/07/2020 10:04

Sample Number: 4
First Solvent: A **Solvent Help:** Grad
Second Solvent: B

Calculation of Gradient accuracy and -precision:

Observed Values for Pumps with Standard Mixing Chamber:

Name	Signal Start mAU	Signal Step mAU	Signal Step 50 mAU	Signal Step 99 mAU	Signal Step 100 mAU
	UV_VIS_1	UV_VIS_1	UV_VIS_1	UV_VIS_1	UV_VIS_1
Solvent change	n.a.	n.a.	n.a.	n.a.	n.a.
STD Gradient_1	0,059	3,466	109,936	215,004	216,248
STD Gradient_2	0,053	3,206	109,743	214,678	215,944
STD Gradient_3	-0,110	3,057	109,494	214,482	215,770

Calculated Steps [%]

	STD A-B	LON	STD C-D	C-D (1.0 ml/min)	Current
Start	0,00	0,00	0,00	0,00	0,00
Step 1	1,00	1,00	0,00	0,00	1,00
Step 2	50,00	99,00	0,00	0,00	50,00
Step 3	99,00	0,00	0,00	0,00	99,00
End	100,00	0,00	0,00	0,00	100,00

Calculation of Ripple:

Observed Values for Pumps: with Standard Mixing Chamber: with Mixing Kit 1 or 2:

Name	Ripple Step 1 mAU	Ripple Step 50 mAU	Ripple Step 99 mAU	Ripple Step 1 mAU	Ripple Step 50 mAU
	UV_VIS_1	UV_VIS_1	UV_VIS_1	UV_VIS_1	UV_VIS_1
Solvent change	n.a.	n.a.	n.a.	n.a.	n.a.
STD Gradient_1	0,090	0,168	0,071	0,157	0,131
STD Gradient_2	0,088	0,095	0,070	0,131	0,084
STD Gradient_3	0,094	0,075	0,067	0,079	0,124

Determination of Pump Unit for Dionex DGPs

Sequence name: PQ_STD_GRAD
 Right end of the sequence GRAD
 Pump's model number: LPG-3400RS
 Pump's model variant: LPG
 Pump unit:

Determination of Pump Flow (Full Qualified and Non-Full-Qualified Variable Name)

CM-Version:	Flow:	
CM6:	2,000	(Pump)
CM7:	2,000	(Pump - DDK driver)
	n.a.	(MicroPump - DDK driver)
	n.a.	(LoadingPump - DDK driver)
	n.a.	(PumpLeft - DDK driver)
	n.a.	(PumpRight - DDK driver)
	2,000	(Pump)
	n.a.	(MicroPump)
	n.a.	(LoadingPump)
	n.a.	(PumpLeft)
	n.a.	(PumpRight)
Used Flow Rate	2,000	(= maximum flow rate)

Determination of Gradient Composition (Full Qualified and Non-Full-Qualified Variab

CM-Version:	%B(0.1):	%B(1.1):	%B(4.1):	%B(8.1):	%B(12.1):
CM6:	n.a.	n.a.	n.a.	n.a.	n.a.
CM7:	0,0	1,0	50,0	99,0	100,0
	n.a.	n.a.	n.a.	n.a.	n.a.
	n.a.	n.a.	n.a.	n.a.	n.a.
	n.a.	n.a.	n.a.	n.a.	n.a.
	n.a.	n.a.	n.a.	n.a.	n.a.
	n.a.	n.a.	n.a.	n.a.	n.a.
Used %B	0,000	1,000	50,000	99,000	100,000

Determination of Gradient Composition (Full Qualified and Non-Full-Qualified Variab

CM-Version:	%B(0.1):	%B(2.1):	%B(8.1):	%B(23.1):	%B(36.1):
CM6:	n.a.	n.a.	n.a.	n.a.	n.a.
CM7:	0,0	1,0	99,0	0,0	0,0
	n.a.	n.a.	n.a.	n.a.	n.a.
	n.a.	n.a.	n.a.	n.a.	n.a.
	n.a.	n.a.	n.a.	n.a.	n.a.
	n.a.	n.a.	n.a.	n.a.	n.a.
	n.a.	n.a.	n.a.	n.a.	n.a.
Used %B	0,000	1,000	99,000	0,000	0,000

Determination of Gradient Composition (Full Qualified and Non-Full-Qualified Variab

CM-Version:	%D(0.1):	%D(1.1):	%D(4.1):	%D(8.1):	%D(12.1):
CM6:	0,0	0,0	0,0	0,0	0,0
CM7:	0,0	0,0	0,0	0,0	0,0
	n.a.	n.a.	n.a.	n.a.	n.a.
	n.a.	n.a.	n.a.	n.a.	n.a.
	n.a.	n.a.	n.a.	n.a.	n.a.
	n.a.	n.a.	n.a.	n.a.	n.a.
	0,0	0,0	0,0	0,0	0,0
Used %D	0,000	0,000	0,000	0,000	0,000

Determination of Gradient Composition (Full Qualified and Non-Full-Qualified Variab

CM-Version:	%D(0.1):	%D(2.1):	%D(8.1):	%D(23.1):	%D(36.1):	
CM6:	0,0	0,0		0,0	0,0	0,0
CM7:	0,0	0,0		0,0	0,0	0,0
	n.a.	n.a.	n.a.	n.a.	n.a.	
	n.a.	n.a.	n.a.	n.a.	n.a.	
	n.a.	n.a.	n.a.	n.a.	n.a.	
	n.a.	n.a.	n.a.	n.a.	n.a.	
	0,0	0,0		0,0	0,0	0,0
Used %D	0,000	0,000	0,000	0,000	0,000	



Gradient_accuracy_8043253 / PUMP_GRADIENT
Printed: 01/07/2020 10:04

Observed Values for Pumps with Mixing Kit 1 or 2:

Signal Start mAU	Signal Step 1 mAU	Signal Step 50 mAU	Signal Step 99 mAU	Signal Step 100 mAU	Signal Step End mAU
UV_VIS_1	UV_VIS_1	UV_VIS_1	UV_VIS_1	UV_VIS_1	UV_VIS_1
n.a.	n.a.	n.a.	n.a.	n.a.	n.a.
0	3,207	109,949	-0,111	n.a.	n.a.
0	3,101	109,752	-0,131	n.a.	n.a.
0	2,978	109,496	-0,239	n.a.	n.a.

with Micro Mixing Kit:

Ripple Step 1 mAU	Ripple Step 50 mAU	Ripple Step 99 mAU
UV_VIS_1	UV_VIS_1	UV_VIS_1
n.a.	n.a.	n.a.
n.a.	0,090	0,168
n.a.	0,088	0,095
n.a.	0,094	0,075
		0,071
		0,070
		0,067

File Name) - STD Gradient A-B

(Pump)
(Pump)
(MicroPump)
(LoadingPump)
(PumpLeft)
(PumpRight)
(Agilent ICF)
(= maximum %B)

File Name) - LON Gradient A-B

(Pump)
(Pump)
(MicroPump)
(LoadingPump)
(PumpLeft)
(PumpRight)
(Agilent ICF)
(= maximum %B)

File Name) - STD Gradient C-D

(Pump)
(Pump)
(MicroPump)
(LoadingPump)
(PumpLeft)
(PumpRight)
(Agilent ICF)
(= maximum %D)

File Name) - Gradient C-D (Flow rate: 1.000 ml/min)

(Pump)
(Pump)
(MicroPump)
(LoadingPump)
(PumpLeft)
(PumpRight)
(Agilent ICF)
(= maximum %B)

Signal Step mAU	Signal Step 1 mAU	Signal Step 50 mAU	Signal Step 99 mAU	Signal Step 100 mAU	Signal End mAU
UV_VIS_1	UV_VIS_1	UV_VIS_1	UV_VIS_1	UV_VIS_1	UV_VIS_1
n.a.	n.a.	n.a.	n.a.	n.a.	n.a.
0,059	3,466	109,936	215,004	216,248	-0,111
0,053	3,206	109,743	214,678	215,944	-0,139
-0,110	3,057	109,494	214,482	215,770	-0,256

Linear Drift correction: $Y = m \cdot X + t$ Drift corrected Signal values: $Y(\text{corr}) = Y(\text{obs}) - m \cdot X - t$

Slope m mAU / min $m = (Y2-Y1)/(X2-X1)$	Offset t mAU $t = Y1 - m \cdot X1$	Signal Step S mAU UV_VIS_1	Signal Step 1 mAU UV_VIS_1	Signal Step 50 mAU UV_VIS_1	Signal Step 99 mAU UV_VIS_1
#VALEUR!	#VALEUR!	#VALEUR!	#VALEUR!	#VALEUR!	#VALEUR!
-0,007499782	0,06682306	0,000	3,429	109,929	215,027
-0,008446109	0,061605985	0,000	3,178	109,749	214,718
-0,006420933	-0,103933589	0,000	3,187	109,650	214,663

Calculated Steps [%]

Signal Step 100 mAU	Signal Step 5 %	Signal Step 1 %	Signal Step 50 %	Signal Step 99 %	Signal Step 100 %
UV_VIS_1	UV_VIS_1	UV_VIS_1	UV_VIS_1	UV_VIS_1	UV_VIS_1
#VALEUR!	#VALEUR!	#VALEUR!	#VALEUR!	#VALEUR!	#VALEUR!
216,293	0,000	1,585	50,824	99,414	100,000
216,009	0,000	1,471	50,808	99,402	100,000
215,970	0,000	1,476	50,771	99,395	100,000



Performance Qualification

• Precision of the Stand. Gradient

• Instruments and Fluidics

Instrument Name	Model	Supplier's Name	Serial Number
Pump	LPG-3400RS	Thermo Scientific	8043253
UV Detector	DAD-3000RS	Thermo Scientific	8043573
Chromeleon Datasystem	V. 6.80 SR16 Build 5387	Thermo Scientific	40640

Accessories	Description
Back Pressure Device	Capillary (L:15 m; ID:0,18 mm)
Solvent A	Water (HPLC-Grade)
Solvent B for Gradient	Water + 0.1% Acetone

• Additional Information

Customer: Customer's Name
 Operator: Operator's Name
 Operator's Jobtitle

Execution Date: 2020-avr-15
 Next Qualification: 2020-oct


• Limits and Test Results

	Limit	Observed max. SD	Result of all Steps
Gradient Precision	<= 0,500 %	0,065 %	Test passed

• Signatures

	Signature	Date
Submitter / Operator:	{seq.submitOperator.userName}	{seq.submitTime}
Reviewer:	{seq.reviewOperator.userName}	{seq.reviewTime}
Approver (e-sig. only):	{seq.approveOperator.userName}	{seq.approveTime}

• Overlay of three Gradients



• Data of three Gradients

Expected Value [%]	Calculated Value Gradient_1 [%]	Calculated Value Gradient_2 [%]	Calculated Value Gradient_3 [%]	Critical SD [%]	Calculated SD [%]	Result
0,00	0,00	0,00	0,00	0,50	0,000	ok
1,00	1,59	1,47	1,48	0,50	0,065	ok
50,00	50,82	50,81	50,77	0,50	0,027	ok
99,00	99,41	99,40	99,39	0,50	0,010	ok
100,00	100,00	100,00	100,00	0,50	0,000	ok

Determination of Pump Unit for Dionex DGPs

Sequence name: PQ_STD_GRAD

Right end of the sequence GRAD

Chromeleon (c) 2017 Thermo Fisher Scientific Inc.
Version 6.80 SR16 Build 5387 (291597)

Gradient_accuracy_8043253 / PUMP_GRADIENT_REPRO

Printed: 01/07/2020 10:04

Pump's model number: LPG-3400RS
Pump's model variant: LPG
Pump unit:

Determination Page Title

Gradient type: • Precision of the Stand. Gradient
Solvent C-D (if used)
Pump unit (DGPs only)
All: • Precision of the Stand. Gradient



