

A18

	A	B	C	D
1	Sample	Sample Name	Ret.Time	Area
2	No.		min	mAU*min
3			Caffeine	Caffeine
4			UV VIS 1	UV VIS 1
5	3	Injector and flow reproducibility_1	1.997	117.7012
6	4	Injector and flow reproducibility_2	2.000	117.2906
7	5	Injector and flow reproducibility_3	2.000	117.4653
8	6	Injector and flow reproducibility_4	2.000	117.8304
9	7	Injector and flow reproducibility_5	2.000	117.8971
10	8	Injector and flow reproducibility_6	2.000	117.7352
11	9	Injector and flow reproducibility_7	2.000	117.9576
12	10	Injector and flow reproducibility_8	1.997	117.8580
13	11	Injector and flow reproducibility_9	2.000	117.6086
14	12	Injector and flow reproducibility_10	2.000	118.0824
15	Average:		1.999	117.7426
16	Rel.Std.Dev:		0.070 %	0.202 %